

# SLB<sup>®</sup> IL (i-series) Capillary GC Columns

Selectivity and Inertness




Improved Inertness

Selectivity Options

Applications

The improved inertness for polar analytes was the inspiration for the development of SLB®-IL (i-series) ionic liquid capillary GC columns. They solve a dilemma that has existed for a long time; is it better to optimize for selectivity **or** inertness? With i-series columns, GC users can enjoy selectivity **and** inertness!





## Improved Inertness

Using the polar column test mix for the quality control of polar and highly polar columns allows an assessment of inertness. For example, the % response of 1-octanol (its peak height relative to a curved line connecting the n-alkane markers) can be measured. A greater value indicates a more inert column.

**Figure 1. Polar Test Mix**

column: SLB-IL76i, 30 m × 0.25 mm I.D., 0.20 μm  
 column: SLB-IL76, 30 m × 0.25 mm I.D., 0.20 μm (28913-U)  
 oven: 130 °C  
 inj. temp.: 250 °C  
 detector: FID, 250 °C  
 carrier gas: helium, 1.2 mL/min (28 cm/sec)  
 injection: 1 μL, 100:1 split  
 liner: 4 mm I.D., split type, cup design  
 sample: polar column test mix (47302),  
 9 analytes, each at 500 μg/mL  
 in methylene chloride

1. *n*-Pentadecane
2. *n*-Hexadecane
3. 2-Octanone
4. *n*-Heptadecane
5. 1-Octanol
6. *n*-Octadecane
7. *n*-Eicosane
8. 2,6-Dimethylaniline
9. 2,6-Dimethylphenol





## Industrial Solvents

SLB-IL (i-series) columns offer selectivity options with inertness for applications involving complex mixtures comprised of compounds with varying functionality.

**Figure 2. Industrial Solvents on SLB-IL60i**

column: SLB-IL60i, 30 m × 0.25 mm I.D., 0.20 μm (29832-U)  
 oven: 40 °C (4 min), 8 °C/min to 200 °C (5 min)  
 inj. temp.: 250 °C  
 detector: FID, 250 °C  
 carrier gas: helium, 30 cm/sec  
 injection: 1 μL, 100:1 split  
 liner: 4 mm I.D., split type, cup design  
 sample: industrial solvents, each at 0.2 % (v/v) in pentane

- |                                       |   |
|---------------------------------------|---|
| 1. Hexane                             | 29. 1,2-Dichloroethane                    |
| 2. 1,1-Dichloroethylene               | 30. <i>n</i> -Butyl acetate               |
| 3. Methyl formate                     | 31. 2-Hexanone                            |
| 4. Acetone                            | 32. Isobutanol                            |
| 5. Ethyl formate                      | 33. Nitropropane                          |
| 6. Methyl acetate                     | 34. Isoamyl acetate                       |
| 7. <i>trans</i> -1,2-Dichloroethylene | 35. Ethylbenzene                          |
| 8. Tetrahydrofuran                    | 36. Mesityl oxide                         |
| 9. Carbon tetrachloride               | 37. <i>p</i> -Xylene                      |
| 10. 1,1-Dichloroethane                | 38. <i>m</i> -Xylene                      |
| 11. Ethyl acetate                     | 39. 5-Methyl-2-hexanone                   |
| 12. Methanol                          | 40. <i>n</i> -Butanol                     |
| 13. Isopropyl acetate                 | 41. <i>n</i> -Amyl acetate                |
| 14. 2-Butanone                        | 42. <i>o</i> -Xylene                      |
| 15. 2-Propanol                        | 43. Isoamyl alcohol                       |
| 16. Methylene chloride                | 44. Chlorobenzene                         |
| 17. Ethanol                           | 45. Styrene                               |
| 18. Benzene                           | 46. 1,1,1,2-Tetrachloroethane             |
| 19. <i>n</i> -Propyl acetate          | 47. Dimethylformamide                     |
| 20. Trichloroethylene                 | 48. Diacetone alcohol                     |
| 21. 4-Methyl-2-pentanone              | 49. Cyclohexanol                          |
| 22. Isobutyl acetate                  | 50. 2-Butoxyethanol<br>(Butyl cellosolve) |
| 23. Tetrachloroethene                 | 51. 1,4-Dichlorobenzene                   |
| 24. Chloroform                        | 52. 1,1,1,2-Tetrachloroethane             |
| 25. <i>sec</i> -Butanol               | 53. 2-Methylphenol                        |
| 26. Toluene                           | 54. 3-Methylphenol                        |
| 27. <i>n</i> -Propanol                | 55. 4-Methylphenol                        |
| 28. 1,4-Dioxane                       |   |





**Figure 3. Industrial Solvents on SLB-IL111i**

column: SLB-IL111i, 30 m × 0.25 mm I.D., 0.20 μm (29883-U)

oven: 40 °C (8 min), 8 °C/min to 200 °C (1 min)

All other conditions and peak IDs are the same as Figure 2.



## Selectivity Options



It is advantageous to have columns with alternative selectivities at hand, because resolution is greatly affected by selectivity. A range of i-series columns were developed and classified as polar (SLB-IL60i), highly polar (SLB-IL76i), and extremely polar (SLB-IL111i). **Table 1** contains complete specifications for all three chemistries.

Two by-products of the fermentation process are active amyl alcohol and isoamyl alcohol. Alcoholic beverage manufacturers are interested in these analytes due to their contribution to the aroma characteristic. The inclusion of n-alkane markers allows some of the selectivity characteristic of each chemistry to be exposed.



**Figure 4. Resolution of Active Amyl Alcohol and Isoamyl Alcohol**

column: SLB-IL60i, 30 m x 0.25 mm I.D., 0.20 μm (29832-U)	detector: FID, 250 °C
column: SLB-IL76i, 30 m x 0.25 mm I.D., 0.20 μm	carrier gas: hydrogen, 40 cm/sec
column: SLB-IL111i, 30 m x 0.25 mm I.D., 0.20 μm (29883-U)	injection: 1 μL, 100:1 split
oven: 90 °C	liner: 4 mm I.D., split type, cup design
inj. temp.: 250 °C	sample: active amyl alcohol, isoamyl alcohol, and n-alkane markers in methylene chloride



### SLB-IL60i



### SLB-IL76i



### SLB-IL111i



SLB-IL (i-series) capillary GC columns are more inert versions of popular ionic liquid chemistries, which provide both selectivity and inertness towards polar analytes, resulting in better accuracy and sensitivity.

**Table 1. SLB-IL (i-series) Column Specifications**

**SLB-IL60i**

**Application:** The selectivity of SLB-IL60i is more polar than PEG/wax phases, resulting in unique elution patterns. It has a higher maximum temperature than most PEG/wax columns (280 °C compared to 260-270 °C). Excellent alternative to existing PEG/wax columns. Also a good GCxGC column choice.

**USP Code:** None

**Phase:** Non-bonded; 1,12-Di(triisopropylphosphonium)dodecane bis(trifluoromethanesulfonyl)imide

**Temp. Limits:** 35 °C to 280 °C (isothermal or programmed)

**SLB-IL76i**

**Application:** The SLB-IL76i phase structure is engineered with numerous interaction mechanisms, resulting in selectivity differences even when compared to columns with similar GC column polarity scale values. Also a good GCxGC column choice.

**USP Code:** None

**Phase:** Non-bonded; Tri(triisopropylphosphoniumhexanamido)triethylamine bis(trifluoromethanesulfonyl)imide

**Temp. Limits:** subambient to 270 °C (isothermal or programmed)

**SLB-IL111i**

**Application:** The selectivity of SLB-IL111i is most orthogonal to non-polar and intermediate polar phases, resulting in very unique elution patterns. Maximum temperature of 270 °C is very impressive for such an extremely polar column. Great choice for separation of polarizable analytes (contain double and/or triple C-C bonds) from neutral analytes. Also a good GCxGC column choice.

**USP Code:** None

**Phase:** Non-bonded; 1,5-Di(2,3-dimethylimidazolium)pentane bis(trifluoromethanesulfonyl)imide

**Temp. Limits:** 50 °C to 260 °C (isothermal or programmed)

**Did you know . . .**

GC stationary phases based on dicationic and polycationic ionic liquids were invented by Prof. Daniel W. Armstrong, currently at the University of Texas at Arlington (USA).

**Ordering Information**

Description	Cat. No.
<b>SLB®-IL60i Capillary GC Columns</b>	
20 m × 0.18 mm I.D., 0.14 μm	<a href="#">29829-U</a>
30 m × 0.25 mm I.D., 0.20 μm	<a href="#">29832-U</a>
60 m × 0.25 mm I.D., 0.20 μm	<a href="#">29833-U</a>
30 m × 0.32 mm I.D., 0.26 μm	<a href="#">29836-U</a>
60 m × 0.32 mm I.D., 0.26 μm	<a href="#">29837-U</a>
<b>SLB®-IL76i Capillary GC Columns</b>	
30 m × 0.25 mm I.D., 0.20 μm	inquire
<b>SLB®-IL111i Capillary GC Columns</b>	
30 m × 0.25 mm I.D., 0.20 μm	<a href="#">29883-U</a>
60 m × 0.25 mm I.D., 0.20 μm	<a href="#">29884-U</a>

**Related Information**

Additional chromatograms, product information, real-time availability, and ordering information is available 24 hours a day at [sigma-aldrich.com/il-gc-inert](http://sigma-aldrich.com/il-gc-inert)

**Order/Customer Service:** [sigma-aldrich.com/order](http://sigma-aldrich.com/order)  
**Technical Service:** [sigma-aldrich.com/techservice](http://sigma-aldrich.com/techservice)  
**Development/Custom Manufacturing Inquiries** **SAFC®** [safcglobal@sial.com](mailto:safcglobal@sial.com)  
**Safety-related Information:** [sigma-aldrich.com/safetycenter](http://sigma-aldrich.com/safetycenter)

3050 Spruce St.  
 St. Louis, MO 63103  
 (314) 771-5765  
[sigma-aldrich.com](http://sigma-aldrich.com)

©2016 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA, SAFC, SIGMA-ALDRICH, ALDRICH and SUPELCO are trademarks of Sigma-Aldrich Co. LLC, registered in the US and other countries. SLB is a registered trademark of Sigma-Aldrich Co. LLC. Sigma-Aldrich, Sigma, Aldrich, Supelco, and SAFC brand products are sold by affiliated Sigma-Aldrich distributors. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at [www.sigmaaldrich.com](http://www.sigmaaldrich.com) and/or on the reverse side of the invoice or packing slip. Sigma-Aldrich Corp. is a subsidiary of Merck KGaA, Darmstadt, Germany.