

Analyses of Underivatized Sugars and Oligosaccharides by HPLC

We tested the typical method of analyzing sugars and oligosaccharides, using HPLC, refractive index detection, and buffered mobile phase. However, for fast, selective analyses we found that HPLC with ultraviolet detection and mobile phases that were not buffered provided shorter analysis times and more accurate separations. The SUPELCOTM LC-NH₂ column is particularly useful for separating these compounds.

Key Words:

- SUPELCOTM columns • sugars • oligosaccharides
- refractive index detection • ultraviolet detection

Although reversed phase HPLC columns based on n-alkyl bonded groups are extremely versatile, they are not suited to every type of analysis. The introduction of functional groups such as nitrile, amino, or diol at the end of the alkyl groups, however, produces new phases which also have useful properties. A SUPELCOTM LC-NH₂ column, packed with a γ -amino-propylsilyl polar bonded phase, is particularly useful for separating sugars and oligosaccharides.


Underivatized monosaccharides are quickly resolved on SUPELCOTM LC-NH₂ columns through normal phase partitioning, based on the interaction of sugar hydroxyl groups with the amino phase (1). Pentoses (i.e., arabinose, xylose) elute before hexoses (i.e., fructose, galactose, glucose, mannose) because pentoses have four hydroxyl groups while hexoses have five (Figure A). The exception is the hexose rhamnose, which elutes with the pentoses because it has only four hydroxyl groups.

For ease in identification, a refractive index (RI) detector is normally used in sugar analyses. However, RI detection generally has limited sensitivity, while temperature variations and pump pulsations must be carefully controlled. In contrast, ultraviolet (UV) detectors are affected much less by temperature and pressure variations. UV detection also enables analysts to use gradient elution to shorten analysis time and separate complex mixtures, including those containing di- and oligosaccharides. Carbohydrate detection is most sensitive at UV wavelengths of 190nm and lower, although solvent impurities can create noise and limit sensitivity at these wavelengths. Figure A compares the use of UV and RI detectors. UV appears to provide the highest sensitivity for D-arabinose, D-fructose, and D-galactose, while D-xylose and α -D-glucose detection is enhanced using RI.

An 85:15 acetonitrile:water mobile phase resolves the most common monosaccharides well (Figure A). Since water and the sugars compete to form hydrogen bonds with the amino phase, sugar retention decreases as water in the mobile phase is increased.

Figure A. Monosaccharides


Column: SUPELCOTM LC-NH₂ (58338), 25cm x 4.6mm ID, 5 μ m particles, with 2cm guard column
Mobile Phase: acetonitrile:water (85:15)
Flow Rate: (A1) 2mL/min
(A2) 1mL/min
Det.: (A1) 190nm UV
(A2) RI
Sample: 10 μ L, component concentrations on figure



Therefore, water content can be increased to shorten analysis time for less complex samples and, similarly, analysts can separate di- and oligosaccharides on SUPELCOTM LC-NH₂ columns by varying the amount of water in the mobile phase. Sorbitol and the common mono- and disaccharides of certain fruit and dietetic products are resolved well on SUPELCOTM LC-NH₂ columns (Figure B). The column also resolves oligosaccharides consisting of 1 to 12 glucose units (DP-1 to DP-12 \bullet). Oligosaccharides consisting of 3 and 4 glucose units usually are resolved poorly on other types of HPLC columns, but they were well resolved in beer samples analyzed on SUPELCOTM LC-NH₂ (Figure C).

Figure B. Food Sugars and Sorbitol

Column: **SUPELCO SIL LC-NH₂, 25cm x 4.6mm ID, 5µm particles**
 Cat. No.: **58338**
 Mobile Phase: acetonitrile:ethanol:water (75:8:17) (v/v/v)
 Flow Rate: 2mL/min
 Temp: 30°C
 Det.: refractive index
 Sample: 20µL (5% w/v each component)




1. Fructose
2. Sorbitol
3. Glucose
4. Sucrose
5. Maltose

797-0603

Figure C. Low Molecular Weight Oligosaccharides in Beer

Column: **SUPELCO SIL LC-NH₂, 25cm x 4.6mm ID, 5µm particles**
 Cat. No.: **58338**
 Mobile Phase: acetonitrile:water (60:40) (v/v)
 Flow Rate: 1mL/min
 Det.: refractive index
 Sample: 10µL of beer



1. Glucose
2. Maltose
3. Maltotriose (DP-3)
4. DP-4
5. DP-5

797-0604

In our work with SUPELCO SIL LC-NH₂ columns, we found that buffered mobile phases, which had been suggested for improving monosaccharide peak symmetry (2), in fact did not appreciably improve peak symmetry. Mobile phases buffered to pH 5-6 generally lengthened retention times without improving resolution. In subsequent analyses without buffers, monosaccharide retention times were shorter than values obtained before buffers were used. When boric acid was used to buffer the mobile phase to pH 5, sugars adsorbed irreversibly to the packing. Consequently, we do not recommend buffered mobile phases in monosaccharide analyses.

SUPELCO SIL LC-NH₂ columns contain a five micron spherical packing. In addition to analytical and guard columns, semi-preparative columns filled with this packing are available. For fast, selective LC analyses of mono-, di-, and oligosaccharides, we highly recommend these aminopropyl-bonded phase columns.

Ordering Information:

SUPELCO SIL LC-NH₂ HPLC Columns

Supplied with fittings to connect to 1/16 inch OD stainless steel tubing.

7.5cm x 4.6mm ID, 3µm particles	58988
15cm x 4.6mm ID, 3µm particles	58989
25cm x 4.6mm ID, 5µm particles	58338

Supelguard LC-NH₂ Guard Column,

2cm x 4.6mm ID, 5µm SUPELCO SIL LC-NH₂ packing

Kit, includes one cartridge, holder, connecting hardware **59558**

Replacement Cartridges, pk. of 2 **59568**

Pelliguard LC-NH₂ Guard Column,

2cm x 4.6mm ID, 40µm pellicular packing

Kit, includes one cartridge, holder, connecting hardware **59646**

Replacement Cartridges, pk. of 4 **59656**

* Degree of polymerization

Pelliguard, SUPELCO SIL, and Supelguard are trademarks of Sigma-Aldrich Co.

References

1. Binder, H., *J. Chromatography*, **189**:414-420 (1980).
2. Rabel, F., et al., *J. Chromatography*, **126**:731-740 (1976).

Acknowledgment

We wish to thank Dr. Eric Coles of the Harvard University Medical School, Boston, for Figure C.

Note 126

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA · Sigma-Aldrich de Argentina, S.A. · Buenos Aires 1119 AUSTRALIA · Sigma-Aldrich Pty. Ltd. · Castle Hill NSW 2154 AUSTRIA · Sigma-Aldrich Handels GmbH · A-1110 Wien
 BELGIUM · Sigma-Aldrich N.V./S.A. · B-2880 Bornem BRAZIL · Sigma-Aldrich Quimica Brasil Ltda. · 01239-010 São Paulo, SP CANADA · Sigma-Aldrich Canada, Ltd. · 2149 Winston Park Dr., Oakville, ON L6H 6J8
 CZECH REPUBLIC · Sigma-Aldrich s.r.o. · 186 00 Praha 8 DENMARK · Sigma-Aldrich Denmark A/S · DK-2665 Vallensbaek Strand FINLAND · Sigma-Aldrich Finland/YA-Kemia Oy · FIN-00700 Helsinki
 FRANCE · Sigma-Aldrich Chimie · 38297 Saint-Quentin-Fallavier Cedex GERMANY · Sigma-Aldrich Chemie GmbH · D-82041 Deisenhofen GREECE · Sigma-Aldrich (o.m.) Ltd. · Ilioupoli 16346, Athens
 HUNGARY · Sigma-Aldrich Kft. · H-1067 Budapest INDIA · Sigma-Aldrich Co. · Bangalore 560 048 IRELAND · Sigma-Aldrich Ireland Ltd. · Dublin 24 ISRAEL · Sigma Israel Chemicals Ltd. · Rehovot 76100
 ITALY · Sigma-Aldrich s.r.l. · 20151 Milano JAPAN · Sigma-Aldrich Japan K.K. · Chuo-ku, Tokyo 103 KOREA · Sigma-Aldrich Korea · Seoul MALAYSIA · Sigma-Aldrich (M) Sdn. Bhd. · Selangor
 MEXICO · Sigma-Aldrich Química S.A. de C.V. · 50200 Toluca NETHERLANDS · Sigma-Aldrich Chemie BV · 3330 AA Zwijndrecht NORWAY · Sigma-Aldrich Norway · Torshov · N-0401 Oslo
 POLAND · Sigma-Aldrich Sp. z o.o. · 61-663 Poznań PORTUGAL · Sigma-Aldrich Química, S.A. · Sintra 2710 RUSSIA · Sigma-Aldrich Russia · Moscow 103062 SINGAPORE · Sigma-Aldrich Pte. Ltd.
 SOUTH AFRICA · Sigma-Aldrich (pty) Ltd. · Jet Park 1459 SPAIN · Sigma-Aldrich Química, S.A. · 28100 Alcobendas, Madrid SWEDEN · Sigma-Aldrich Sweden AB · 135 70 Stockholm
 SWITZERLAND · Supelco · CH-9471 Buchs UNITED KINGDOM · Sigma-Aldrich Company Ltd. · Poole, Dorset BH12 4QH
 UNITED STATES · Supelco · Supelco Park · Bellefonte, PA 16823-0048 · Phone 800-247-6628 or 814-359-3441 · Fax 800-447-3044 or 814-359-3044 · email: supelco@sial.com

H

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

BIA